

Super-spherical markers

High quality threaded markers.



Product information

Qualisys super-spherical markers are available in 10 different variants. The markers are designed for maximum roundness and with their solid build, they are more durable than their hollow counterparts.

The markers are available in packs of 10, 50 and 100 pieces, including plastic bases. Prices below do not include taxes and shipping costs which vary by location.

To get a price including these costs, please contact your local distributor or sales@qualisys.com.

Replacement plastic bases are also available separately.



Underwater markers



	19 mm	16 mm	14 mm	12.5 mm	9.5 mm	8 mm	6.5 mm	19 mm	16 mm	14 mm
Size	19 mm	16 mm	14 mm	12.5 mm	9.5 mm	8 mm	6.5 mm	19 mm	16 mm	14 mm
Weight	2.5 g	2.0 g	1.5 g	1.0 g	0.5 g	0.3 g	0.2 g	2.5 g	2.0 g	1.5 g
Thread	M4	M4	M4	M4	M3	M3	M3	M4	M4	M4
10 pack	€55	€55	€55	€55	€55	€75	€75	€55	€55	€55
50 pack	€240	€240	€240	€240	€240	€335	€335	€240	€240	€240
100 pack	€445	€445	€445	€445	€445	€610	€610	€445	€445	€445

Read about our other markers on www.qualisys.com/products/accessories/.